

产品名称: CD272 Rabbit Monoclonal Antibody
产品货号: AMRe01784



产品概述 (Summary)

产品名称 (Production Name)	CD272 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Btla
别名 (Alternative Names)	A630002H24; B-and T-lymphocyte-associated protein
基因 ID (Gene ID)	208154.0
蛋白 ID (SwissProt ID)	Q7TSA3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 55-80 kDa

研究背景 (Background)

Lymphocyte inhibitory receptor which inhibits lymphocytes during immune response.


产品名称: CD272 Rabbit Monoclonal Antibody
产品货号: AMRe01784




研究领域 (Research Area)

Immunology


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded mouse spleen using CD272 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD272/BTLA in 3T3 lysates using CD272/BTLA antibody.



Western blot analysis of CD272/BTLA in K562, HeLa lysates using CD272/BTLA antibody

产品名称: CD272 Rabbit Monoclonal Antibody
产品货号: AMRe01784



注意事项 (Note)

For research use only .