

产品名称: CD276 Rabbit Monoclonal Antibody
产品货号: AMRe01785



产品概述 (Summary)

产品名称 (Production Name)	CD276 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cd276
别名 (Alternative Names)	B7h3; B7RP-2; AU016588; 6030411F23Rik
基因 ID (Gene ID)	102657.0
蛋白 ID (SwissProt ID)	Q8VE98.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 90-110 kDa

研究背景 (Background)

Modulates T-cell-mediated immune responses and the development of acute and chronic transplant rejection. Plays a

产品名称: CD276 Rabbit Monoclonal Antibody
产品货号: AMRe01785




positive regulatory role in bone formation and has a dual role in the bone-immune interface. Induces antitumor immunity as it activates both acquired and innate immunity leading to natural killer cell and CD8 T-cell dependent killing of tumor cells.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD276 in HepG2 lysates using CD276 antibody.

注意事项 (Note)

For research use only .