

产品名称: Phospho-CD3 zeta (Tyr142) Rabbit Monoclonal Antibody
产品货号: AMRe01788



产品概述 (Summary)

产品名称 (Production Name)	Phospho-CD3 zeta (Tyr142) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD247
别名 (Alternative Names)	CD247; CD3Z; T3Z; TCRZ; T-cell surface glycoprotein CD3 zeta chain; T-cell receptor T3 zeta chain; CD antigen CD247
基因 ID (Gene ID)	919
蛋白 ID (SwissProt ID)	P20963.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 19 kDa; Observed MW: 19 kDa

产品名称: Phospho-CD3 zeta (Tyr142) Rabbit Monoclonal Antibody
产品货号: AMRe01788

研究背景 (Background)


Defects in CD3D are a cause of severe combined immunodeficiency autosomal recessive

T-cell-negative/B-cell-positive/NK-cell-positive (T(-)/B(+)/NK(+)) SCID [MIM:608971]. A form of severe combined immunodeficiency (SCID), a genetically and clinically heterogeneous group of rare congenital disorders characterized by impairment of both humoral and cell-mediated immunity, leukopenia, and low or absent antibody levels.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of Phospho-CD3 zeta (Tyr142) in Jurkat lysates using Phospho-CD3 zeta (Tyr142) antibody.

注意事项 (Note)

For research use only .