

产品名称: CD35 Rabbit Monoclonal Antibody
产品货号: AMRe01791



产品概述 (Summary)

产品名称 (Production Name)	CD35 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CR1
别名 (Alternative Names)	C3 binding protein; C3b/C4b receptor; C3BR; C4BR; CD35; CR1; KN
基因 ID (Gene ID)	1378
蛋白 ID (SwissProt ID)	P17927.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 224 kDa; Observed MW: 280 kDa

研究背景 (Background)

Mediates cellular binding of particles and immune complexes that have activated complement.


产品名称: CD35 Rabbit Monoclonal Antibody
产品货号: AMRe01791



研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of CD35 in Human tonsil lysates using CD35 antibody.

注意事项 (Note)

For research use only .