

产品名称: CD39 Rabbit Monoclonal Antibody
产品货号: AMRe01792



产品概述 (Summary)

产品名称 (Production Name)	CD39 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ENTPD1 ENTPD1; CD39; Ectonucleoside triphosphate diphosphohydrolase 1; NTPDase 1; Ecto-ATP diphosphohydrolase 1; Ecto-ATPDase 1; Ecto-ATPase 1; Ecto-apyrase; Lymphoid cell activation antigen; CD antigen CD39
别名 (Alternative Names)	
基因 ID (Gene ID)	953
蛋白 ID (SwissProt ID)	P49961.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 58 kDa; Observed MW: 34-70 kDa

产品名称: CD39 Rabbit Monoclonal Antibody
产品货号: AMRe01792




研究背景 (Background)

In the nervous system, could hydrolyze ATP and other nucleotides to regulate purinergic neurotransmission. Could also be implicated in the prevention of platelet aggregation. Hydrolyzes ATP and ADP equally well.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD39 in mouse spleen lysates using CD39 antibody.

注意事项 (Note)

For research use only .