

产品名称: CD99 Rabbit Monoclonal Antibody
产品货号: AMRe01801



产品概述 (Summary)

产品名称 (Production Name)	CD99 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD99
别名 (Alternative Names)	MIC2; HBA71; MIC2X; MIC2Y; MSK5X
基因 ID (Gene ID)	4267
蛋白 ID (SwissProt ID)	P14209.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 19 kDa; Observed MW: 30 kDa

研究背景 (Background)

Involved in T-cell adhesion processes. It is involved in spontaneous rosette formation with erythrocytes.


产品名称: CD99 Rabbit Monoclonal Antibody
产品货号: AMRe01801




研究领域 (Research Area)

Immunology


图片 (Image Data)



Western blot analysis of CD99 in MOLT4 lysates using CD99 antibody.



Western blot analysis of CD99 in 3T3 lysates using CD99 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using CD99 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: CD99 Rabbit Monoclonal Antibody
产品货号: AMRe01801



注意事项 (Note)

For research use only .