

产品名称: CDK9 Rabbit Monoclonal Antibody
产品货号: AMRe01810



产品概述 (Summary)

产品名称 (Production Name)	CDK9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDK9
别名 (Alternative Names)	TAK; C-2k; CTK1; CDC2L4; PITALRE
基因 ID (Gene ID)	1025
蛋白 ID (SwissProt ID)	P50750.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa; Observed MW: 43 kDa

研究背景 (Background)

Cyclin dependent kinases (CDKs) are activated in part by cyclin binding and by phosphorylation of a conserved threonine in

产品名称: CDK9 Rabbit Monoclonal Antibody
产品货号: AMRe01810




the T-loop domain. Member of the cyclin-dependent kinase pair (CDK9/cyclin-T) complex, also called positive transcription elongation factor b (P-TEFb), which facilitates the transition from abortive to production elongation by phosphorylating the CTD (C-terminal domain) of the large subunit of RNA polymerase II (RNAP II), SUPT5H and RDBP.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of Cdk9 in K562, C6, 3T3, HeLa lysates using CDK9 antibody.



Immunocytochemistry analysis of Cdk9(green) in HeLa using Cdk9 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Cdk9 antibody. High-pressure and

产品名称: CDK9 Rabbit Monoclonal Antibody
产品货号: AMRe01810



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .