

产品名称: CENPE Rabbit Monoclonal Antibody
产品货号: AMRe01818



产品概述 (Summary)

产品名称 (Production Name)	CENPE Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CENPE
别名 (Alternative Names)	CENP E; Centromere associated protein E; Centromere protein E 312kDa; KIF10; Kinesin family member 10; Kinesin related protein; PPP1R61
基因 ID (Gene ID)	1062
蛋白 ID (SwissProt ID)	Q02224.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 316 kDa; Observed MW: 316 kDa

研究背景 (Background)

产品名称: CENPE Rabbit Monoclonal Antibody
产品货号: AMRe01818




CENPE is a 250-300 kDa human centromere-associated kinesin-like motor protein that accumulates in G2 phase.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of CENPE in Jurkat lysates using CENPE antibody.

注意事项 (Note)

For research use only .