

产品名称: CGBP Rabbit Monoclonal Antibody
产品货号: AMRe01820



产品概述 (Summary)

产品名称 (Production Name)	CGBP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CXXC1
别名 (Alternative Names)	CFP1; CGBP; SPP1; PCCX1; PHF18; hCGBP; ZCGPC1; HsT2645; 2410002I16Rik; 5830420C16Rik
基因 ID (Gene ID)	30827
蛋白 ID (SwissProt ID)	Q9P0U4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 76 kDa; Observed MW: 76 kDa

研究背景 (Background)

产品名称: CGBP Rabbit Monoclonal Antibody
产品货号: AMRe01820




Transcriptional activator that exhibits a unique DNA binding specificity for CpG unmethylated motifs with a preference for CpGG.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of CGBP in C6, 3T3, HeLa lysates using CGBP antibody.



Immunocytochemistry analysis of CGBP(green) in HeLa using CGBP antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using CGBP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: CGBP Rabbit Monoclonal Antibody
产品货号: AMRe01820



注意事项 (Note)

For research use only .