

产品名称: cIAP1 Rabbit Monoclonal Antibody
产品货号: AMRe01826



产品概述 (Summary)

产品名称 (Production Name)	cIAP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BIRC2
别名 (Alternative Names)	cIAP1; C-IAP1; IAP homolog B; Inhibitor of apoptosis protein 2; IAP-2; hIAP-2; hIAP2; API1; IAP2; MIHB; RNF48
基因 ID (Gene ID)	329
蛋白 ID (SwissProt ID)	Q13490.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 70 kDa

研究背景 (Background)

产品名称: cIAP1 Rabbit Monoclonal Antibody
产品货号: AMRe01826




Apoptotic suppressor. The BIR motifs region interacts with TNF receptor associated factors 1 and 2 (TRAF1 and TRAF2) to form an heteromeric complex, which is then recruited to the tumor necrosis factor receptor 2 (TNFR2).


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of cIAP1 in K562 lysates using cIAP1 antibody.



Immunocytochemistry analysis of cIAP1(green) in HeLa using cIAP1 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .