

产品名称: CRMP3 Rabbit Monoclonal Antibody
产品货号: AMRe01854



产品概述 (Summary)

产品名称 (Production Name)	CRMP3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DPYSL4 DPYSL4; CRMP3; ULIP4; Dihydropyrimidinase-related protein 4; DRP-4;
别名 (Alternative Names)	Collapsin response mediator protein 3; CRMP-3; UNC33-like phosphoprotein 4; ULIP-4
基因 ID (Gene ID)	10570
蛋白 ID (SwissProt ID)	O14531.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa; Observed MW: 62 kDa

产品名称: CRMP3 Rabbit Monoclonal Antibody
产品货号: AMRe01854




研究背景 (Background)

Necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton. Plays a role in axon guidance, neuronal growth cone collapse and cell migration .

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of CRMP3 in rat Brain lysates using CRMP3 antibody.

注意事项 (Note)

For research use only .