

产品名称: Cytokeratin 19 Rabbit Monoclonal Antibody
产品货号: AMRe01888



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 19 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT19
别名 (Alternative Names)	KRT19; Keratin; type I cytoskeletal 19; Cytokeratin-19; CK-19; Keratin-19; K19
基因 ID (Gene ID)	3880
蛋白 ID (SwissProt ID)	P08727.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa; Observed MW: 41 kDa

研究背景 (Background)

Involved in the organization of myofibers. Together with KRT8, helps to link the contractile apparatus to dystrophin at the

产品名称: Cytokeratin 19 Rabbit Monoclonal Antibody
产品货号: AMRe01888




costameres of striated muscle.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of Cytokeratin 19 in rat Brain lysates using Cytokeratin 19 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Cytokeratin 19 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Cytokeratin 19 (green) in MCF-7 using Cytokeratin 19 antibody, and DAPI (blue)

产品名称: Cytokeratin 19 Rabbit Monoclonal Antibody
产品货号: AMRe01888



注意事项 (Note)

For research use only .