

产品名称: DNAJC8 Rabbit Monoclonal Antibody
产品货号: AMRe01913



产品概述 (Summary)

产品名称 (Production Name)	DNAJC8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNAJC8
别名 (Alternative Names)	SPF31; HSPC331
基因 ID (Gene ID)	22826
蛋白 ID (SwissProt ID)	O75937.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 30 kDa; Observed MW: 30 kDa

研究背景 (Background)

Suppresses polyglutamine (polyQ) aggregation of ATXN3 in neuronal cells (PubMed:27133716).


产品名称: DNAJC8 Rabbit Monoclonal Antibody
产品货号: AMRe01913




研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of DjC8 in mouse brain, rat lung lysates using DNAJC8 antibody.



Western blot analysis of DjC8 in Jurkat, HeLa, COS7, Raji lysates using DjC8 antibody.

注意事项 (Note)

For research use only .