

产品名称: DKK3 Rabbit Monoclonal Antibody
产品货号: AMRe01914



产品概述 (Summary)

产品名称 (Production Name)	DKK3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DKK3
别名 (Alternative Names)	DKK3; REIC; Dickkopf-related protein 3; Dickkopf-3; Dkk-3; hDkk-3
基因 ID (Gene ID)	27122
蛋白 ID (SwissProt ID)	Q9UBP4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa; Observed MW: 70 kDa

研究背景 (Background)

Antagonizes canonical Wnt signaling by inhibiting LRP5/6 interaction with Wnt and by forming a ternary complex with the

产品名称: DKK3 Rabbit Monoclonal Antibody
产品货号: AMRe01914




transmembrane protein KREMEN that promotes internalization of LRP5/6.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Dkk3 in rat Brain lysates using DKK3 antibody.

注意事项 (Note)

For research use only .