

产品名称：DLX4 Rabbit Monoclonal Antibody
产品货号：AMRe01915



产品概述 (Summary)

产品名称 (Production Name)	DLX4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DLX4
别名 (Alternative Names)	DLX4; BP1; DLX7; DLX8; DLX9; Homeobox protein DLX-4; Beta protein 1; Homeobox protein DLX-7; Homeobox protein DLX-8
基因 ID (Gene ID)	1748
蛋白 ID (SwissProt ID)	Q92988.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 35 kDa

研究背景 (Background)

产品名称: DLX4 Rabbit Monoclonal Antibody
产品货号: AMRe01915




May play a role in determining the production of hemoglobin S. May act as a repressor. During embryonic development, plays a role in palatogenesis.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Immunocytochemistry analysis of DLX4 (green) in K562 using DLX4 antibody, and DAPI(blue).



Western blot analysis of Dlx4 in K562, rat Brain, Hela lysates using Dlx4 antibody.

注意事项 (Note)

For research use only .