

产品名称: E1B AP5 Rabbit Monoclonal Antibody
产品货号: AMRe01927



产品概述 (Summary)

产品名称 (Production Name)	E1B AP5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HNRNPUL1
别名 (Alternative Names)	E1BAP5; E1B-AP5; HNRPUL1
基因 ID (Gene ID)	11100
蛋白 ID (SwissProt ID)	Q9BUJ2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 96 kDa; Observed MW: 96 kDa

研究背景 (Background)

Acts as a basic transcriptional regulator. Represses basic transcription driven by several virus and cellular promoters. When

产品名称: E1B AP5 Rabbit Monoclonal Antibody
产品货号: AMRe01927




associated with BRD7, activates transcription of glucocorticoid-responsive promoter in the absence of ligand-stimulation. Plays also a role in mRNA processing and transport. Binds avidly to poly(G) and poly(C) RNA homopolymers in vitro.

研究领域 (Research Area)

Microbiology

图片 (Image Data)



Western blot analysis of E1BAP5 in K562, C6 lysates using E1B AP5 antibody.

注意事项 (Note)

For research use only .