

产品名称: ELMO1 Rabbit Monoclonal Antibody
产品货号: AMRe01944



产品概述 (Summary)

产品名称 (Production Name)	ELMO1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ELMO1
别名 (Alternative Names)	ELMO1; KIAA0281; Engulfment and cell motility protein 1; Protein ced-12 homolog
基因 ID (Gene ID)	9844
蛋白 ID (SwissProt ID)	Q92556.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 84 kDa; Observed MW: 84 kDa

研究背景 (Background)

产品名称: ELMO1 Rabbit Monoclonal Antibody
产品货号: AMRe01944




Involved in cytoskeletal rearrangements required for phagocytosis of apoptotic cells and cell motility. Acts in association with DOCK1 and CRK.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Immunocytochemistry analysis of ELMO1 (green) in Jurkat using ELMO1 antibody, and DAPI (blue).



Western blot analysis of ELMO1 in Jurkat lysates using ELMO1 antibody.



Western blot analysis of ELMO1 in rat Brain, C6, 3T3 lysates using ELMO1 antibody

产品名称: ELMO1 Rabbit Monoclonal Antibody
产品货号: AMRe01944



注意事项 (Note)

For research use only .