

产品名称: Emerin Rabbit Monoclonal Antibody
产品货号: AMRe01946



产品概述 (Summary)

产品名称 (Production Name)	Emerin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EMD
别名 (Alternative Names)	STA; EDMD; LEMD5
基因 ID (Gene ID)	2010
蛋白 ID (SwissProt ID)	P50402.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 35 kDa

研究背景 (Background)

Emerin may regulate gene expression through binding to other transcriptional regulators. Emerin binds to β -catenin and

产品名称: Emerin Rabbit Monoclonal Antibody
产品货号: AMRe01946




inhibits its nuclear accumulation. Recent studies demonstrate that Emerin is required for HIV-1 infectivity.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of Emerin in 3T3, HeLa lysates using Emerin antibody.



Immunocytochemistry analysis of Emerin(green) in HeLa using Emerin antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Emerin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Emerin Rabbit Monoclonal Antibody
产品货号: AMRe01946



注意事项 (Note)

For research use only .