

产品名称: ETAA1 Rabbit Monoclonal Antibody
产品货号: AMRe01961



产品概述 (Summary)

产品名称 (Production Name)	ETAA1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ETAA1
别名 (Alternative Names)	ETAA16
基因 ID (Gene ID)	54465
蛋白 ID (SwissProt ID)	Q9NY74.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 103 kDa; Observed MW: 103 kDa

研究背景 (Background)

Replication stress response protein that accumulates at DNA damage sites and promotes replication fork progression and

产品名称: ETAA1 Rabbit Monoclonal Antibody
产品货号: AMRe01961




integrity (PubMed:27601467, PubMed:27723720, PubMed:27723717). Recruited to stalled replication forks via interaction with the RPA complex and directly stimulates ATR kinase activity independently of TOPBP1 (PubMed:27723720, PubMed:27723717). Probably only regulates a subset of ATR targets (PubMed:27723720, PubMed:27723717).

研究领域 (Research Area)

Cancer

图片 (Image Data)



Western blot analysis of ETAA1 in K562 lysates using ETAA1 antibody.

注意事项 (Note)

For research use only .