

产品名称: EWSR1 Rabbit Monoclonal Antibody
产品货号: AMRe01962



产品概述 (Summary)

产品名称 (Production Name)	EWSR1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EWSR1
别名 (Alternative Names)	EWSR1; EWS; RNA-binding protein EWS; EWS oncogene; Ewing sarcoma breakpoint region 1 protein
基因 ID (Gene ID)	2130
蛋白 ID (SwissProt ID)	Q01844.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa; Observed MW: 85 kDa

研究背景 (Background)

产品名称: EWSR1 Rabbit Monoclonal Antibody
产品货号: AMRe01962




Might normally function as a repressor. EWS-fusion-proteins (EFPS) may play a role in the tumorigenic process. They may disturb gene expression by mimicking, or interfering with the normal function of CTD-POLII within the transcription initiation complex.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Immunocytochemistry analysis of EWSR1 (green) in 293t using EWSR1 antibody, and DAPI(blue).



Western blot analysis of EWSR1 in HeLa lysates using EWSR1 antibody.

注意事项 (Note)

For research use only .