

产品名称: EWSR1 Rabbit Monoclonal Antibody
产品货号: AMRe01963



产品概述 (Summary)

产品名称 (Production Name)	EWSR1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EWSR1
别名 (Alternative Names)	EWSR1; EWS; RNA-binding protein EWS; EWS oncogene; Ewing sarcoma breakpoint region 1 protein
基因 ID (Gene ID)	2130
蛋白 ID (SwissProt ID)	Q01844.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa; Observed MW: 85 kDa

研究背景 (Background)

产品名称: EWSR1 Rabbit Monoclonal Antibody
产品货号: AMRe01963




Might normally function as a repressor. EWS-fusion-proteins (EFPS) may play a role in the tumorigenic process. They may disturb gene expression by mimicking, or interfering with the normal function of CTD-POLII within the transcription initiation complex.


研究领域 (Research Area)

Tags & Cell Markers


图片 (Image Data)



Western blot analysis of EWSR1 in rat brain lysates using EWSR1 antibody.



Western blot analysis of EWSR1 in MCF-7, 3T3, HeLa lysates using EWSR1 antibody.



Immunocytochemistry analysis of EWSR1(green) in HeLa using EWSR1 antibody, and DAPI(blue)

产品名称: EWSR1 Rabbit Monoclonal Antibody
产品货号: AMRe01963



注意事项 (Note)

For research use only .