

产品名称: FAM98B Rabbit Monoclonal Antibody
产品货号: AMRe01973



产品概述 (Summary)

产品名称 (Production Name)	FAM98B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FAM98B
别名 (Alternative Names)	Protein FAM98B
基因 ID (Gene ID)	283742
蛋白 ID (SwissProt ID)	Q52LJ0.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 46 kDa; Observed MW: 37 kDa

研究背景 (Background)

Positively stimulates PRMT1-induced protein arginine dimethylated arginine methylation (PubMed:28040436). Promotes

产品名称: FAM98B Rabbit Monoclonal Antibody
产品货号: AMRe01973




colorectal cancer cell malignancy (PubMed:28040436).


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of FAM98B in K562 lysates using FAM98B antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FAM98B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of FAM98B(green) in HeLa using FAM98B antibody, and DAPI(blue)

产品名称: FAM98B Rabbit Monoclonal Antibody
产品货号: AMRe01973



注意事项 (Note)

For research use only .