

产品名称: FBP1 Rabbit Monoclonal Antibody
产品货号: AMRe01976



产品概述 (Summary)

产品名称 (Production Name)	FBP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FBP1
别名 (Alternative Names)	FBP; FBP1; FBPase 1
基因 ID (Gene ID)	2203
蛋白 ID (SwissProt ID)	P09467.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa; Observed MW: 37 kDa

研究背景 (Background)

Catalyzes the hydrolysis of fructose 1,6-bisphosphate to fructose 6-phosphate in the presence of divalent cations, acting as

产品名称: FBP1 Rabbit Monoclonal Antibody
产品货号: AMRe01976




a rate-limiting enzyme in gluconeogenesis.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of FBP1 in rat Brain lysates using FBP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .