

产品名称: FCGR1A Rabbit Monoclonal Antibody
产品货号: AMRe01978



产品概述 (Summary)

产品名称 (Production Name)	FCGR1A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FCGR1A
别名 (Alternative Names)	Fc-gamma RI; FcRI; Fc-gamma R1A; FcgammaR1a; CD64
基因 ID (Gene ID)	2209
蛋白 ID (SwissProt ID)	P12314.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa; Observed MW: 43-80 kDa

研究背景 (Background)

High affinity receptor for the Fc region of immunoglobulins gamma. Functions in both innate and adaptive immune

产品名称: FCGR1A Rabbit Monoclonal Antibody
产品货号: AMRe01978




responses.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of FCGR1A in THP-1 lysates using FCGR1A antibody.

注意事项 (Note)

For research use only .