

产品名称: FNDC5 Rabbit Monoclonal Antibody
产品货号: AMRe01991



产品概述 (Summary)

产品名称 (Production Name)	FNDC5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FNDC5
别名 (Alternative Names)	FNDC 5; FRCP2; irisin
基因 ID (Gene ID)	252995
蛋白 ID (SwissProt ID)	Q8NAU1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa; Observed MW: 21 kDa

研究背景 (Background)

Irisin: Contrary to mouse, may not be involved in the beneficial effects of muscular exercise, nor in the induction of

产品名称: FNDC5 Rabbit Monoclonal Antibody
产品货号: AMRe01991




browning of human white adipose tissue.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of FNDC5 in K562, C6, 3T3, HeLa lysates using FNDC5 antibody.

注意事项 (Note)

For research use only .