

产品名称: FNTB Rabbit Monoclonal Antibody
产品货号: AMRe01993



产品概述 (Summary)

产品名称 (Production Name)	FNTB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FNTB
别名 (Alternative Names)	FNTB; Protein farnesyltransferase subunit beta; FTase-beta; CAAX farnesyltransferase subunit beta; Ras proteins prenyltransferase subunit beta
基因 ID (Gene ID)	100529261
蛋白 ID (SwissProt ID)	P49356.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 49 kDa; Observed MW: 49 kDa

研究背景 (Background)

产品名称: FNTB Rabbit Monoclonal Antibody
产品货号: AMRe01993




Essential subunit of the farnesyltransferase complex. Catalyzes the transfer of a farnesyl moiety from farnesyl diphosphate to a cysteine at the fourth position from the C-terminus of several proteins having the C-terminal sequence Cys-aliphatic-aliphatic-X.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of FNTB in K562, HeLa lysates using FNTB antibody.

注意事项 (Note)

For research use only .