

产品名称: FOXP1 Rabbit Monoclonal Antibody
产品货号: AMRe02000



产品概述 (Summary)

产品名称 (Production Name)	FOXP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOXP1
别名 (Alternative Names)	FOXP1; HSPC215; Forkhead box protein P1
基因 ID (Gene ID)	27086
蛋白 ID (SwissProt ID)	Q9H334.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 75 kDa; Observed MW: 82-90 kDa

研究背景 (Background)

Transcriptional repressor. It plays an important role in the specification and differentiation of lung epithelium. Can act with

产品名称: FOXP1 Rabbit Monoclonal Antibody
产品货号: AMRe02000




CTBP1 to synergistically repress transcription but CTPBP1 is not essential. Essential transcriptional regulator of B cell development.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of FOXP1 in Raw264.7, rat lung lysates using FOXP1 antibody.

注意事项 (Note)

For research use only .