

产品名称: FRAT1 Rabbit Monoclonal Antibody
产品货号: AMRe02004



产品概述 (Summary)

产品名称 (Production Name)	FRAT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FRAT1
别名 (Alternative Names)	FRAT-1
基因 ID (Gene ID)	10023
蛋白 ID (SwissProt ID)	Q92837.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 29 kDa

研究背景 (Background)

Positively regulates the Wnt signaling pathway by stabilizing beta-catenin through the association with GSK-3. May play a

产品名称: FRAT1 Rabbit Monoclonal Antibody
产品货号: AMRe02004




role in tumor progression and collaborate with PIM1 and MYC in lymphomagenesis.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Immunocytochemistry analysis of FRAT1 (green) in LNcap using FRAT1 antibody, and DAPI(blue).



Western blot analysis of FRAT1 in K562, HeLa lysates using FRAT1 antibody.

注意事项 (Note)

For research use only .