

产品名称: FUBP1 Rabbit Monoclonal Antibody
产品货号: AMRe02009



产品概述 (Summary)

产品名称 (Production Name)	FUBP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FUBP1
别名 (Alternative Names)	DNA helicase V; FBP; FUBP; Fubp1; Fubp4; FUSE-binding protein 1; HDH V
基因 ID (Gene ID)	8880
蛋白 ID (SwissProt ID)	Q96AE4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa; Observed MW: 74 kDa

研究背景 (Background)

Regulates MYC expression by binding to a single-stranded far-upstream element upstream of the MYC promoter. May act

产品名称: FUBP1 Rabbit Monoclonal Antibody
产品货号: AMRe02009




both as activator and repressor of transcription.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Immunocytochemistry analysis of FUBP1 (green) in 293 using FUBP1 antibody, and DAPI (blue).




Western blot analysis of FUBP1/GBP in mouse brain, mouse spleen lysates using FUBP1/GBP antibody.



Western blot analysis of FUBP1/GBP in Jurkat, C6 lysates using FUBP1/GBP antibody

产品名称: FUBP1 Rabbit Monoclonal Antibody
产品货号: AMRe02009



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using FUBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .