

产品名称: GCET2 Rabbit Monoclonal Antibody
产品货号: AMRe02023



产品概述 (Summary)

产品名称 (Production Name)	GCET2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GCSAM
别名 (Alternative Names)	GCAT2; Gcet; Gcsam; hGAL; M17; M17 L
基因 ID (Gene ID)	257144
蛋白 ID (SwissProt ID)	Q8N6F7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 21 kDa

研究背景 (Background)

Involved in the negative regulation of lymphocyte motility. It mediates the migration-inhibitory effects of IL6. Serves as a

产品名称: GCET2 Rabbit Monoclonal Antibody
产品货号: AMRe02023




positive regulator of the RhoA signaling pathway. Enhancement of RhoA activation results in inhibition of lymphocyte and lymphoma cell motility by activation of its downstream effector ROCK. Is a regulator of B-cell receptor signaling, that acts through SYK kinase activation.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of GCET2 in Raji lysates using GCET2 antibody.

注意事项 (Note)

For research use only .