

产品名称: Gelsolin Rabbit Monoclonal Antibody
产品货号: AMRe02027



产品概述 (Summary)

产品名称 (Production Name)	Gelsolin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSN
别名 (Alternative Names)	ADF; AGEL; Brevin; Gelsolin; Gsn
基因 ID (Gene ID)	2934
蛋白 ID (SwissProt ID)	P06396.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 86 kDa; Observed MW: 86 kDa

研究背景 (Background)

Calcium-regulated, actin-modulating protein that binds to the plus (or barbed) ends of actin monomers or filaments,

产品名称: Gelsolin Rabbit Monoclonal Antibody
产品货号: AMRe02027




preventing monomer exchange (end-blocking or capping). It can promote the assembly of monomers into filaments (nucleation) as well as sever filaments already formed. Plays a role in ciliogenesis.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Gelsolin in MCF-7, C6, HeLa lysates using Gelsolin antibody.

注意事项 (Note)

For research use only .