

产品名称: GFPT1 Rabbit Monoclonal Antibody
产品货号: AMRe02032



产品概述 (Summary)

产品名称 (Production Name)	GFPT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GFPT1
别名 (Alternative Names)	GFA; GFAT; GFPT; MSLG; CMS12; GFAT1; CMSTA1; GFAT 1; GFAT1m; GFPT1L
基因 ID (Gene ID)	2673
蛋白 ID (SwissProt ID)	Q06210.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 79 kDa; Observed MW: 79 kDa

研究背景 (Background)

Controls the flux of glucose into the hexosamine pathway. Most likely involved in regulating the availability of precursors

产品名称: GFPT1 Rabbit Monoclonal Antibody
产品货号: AMRe02032




for N- and O-linked glycosylation of proteins. Regulates the circadian expression of clock genes ARNTL/BMAL1 and CRY1.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of GFPT1 in HeLa lysates using GFPT1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using GFPT1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .