

产品名称: GFPT2 Rabbit Monoclonal Antibody
产品货号: AMRe02033



产品概述 (Summary)

产品名称 (Production Name)	GFPT2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GFPT2
别名 (Alternative Names)	GFAT; GFAT2; GFAT 2
基因 ID (Gene ID)	9945
蛋白 ID (SwissProt ID)	O94808.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 77 kDa; Observed MW: 77 kDa

研究背景 (Background)

Controls the flux of glucose into the hexosamine pathway. Most likely involved in regulating the availability of precursors

产品名称: GFPT2 Rabbit Monoclonal Antibody
产品货号: AMRe02033




for N- and O-linked glycosylation of proteins.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Immunocytochemistry analysis of GFPT2 (green) in HT-1080 using GFPT2 antibody, and DAPI(blue).



Western blot analysis of GFPT2 in rat spleen lysates using GFPT2 antibody.



Western blot analysis of GFPT2 in 3T3, HeLa lysates using GFPT2 antibody

产品名称: GFPT2 Rabbit Monoclonal Antibody
产品货号: AMRe02033



注意事项 (Note)

For research use only .