

产品名称: Glucocorticoid Receptor Rabbit Monoclonal Antibody
产品货号: AMRe02037



产品概述 (Summary)

产品名称 (Production Name)	Glucocorticoid Receptor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Nr3C1
别名 (Alternative Names)	GR; GCR; GRL; GCCR; GCRST
基因 ID (Gene ID)	14815.0
蛋白 ID (SwissProt ID)	P06537.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 87 kDa; Observed MW: 94,91 kDa

研究背景 (Background)

产品名称: Glucocorticoid Receptor Rabbit Monoclonal Antibody
产品货号: AMRe02037




Receptor for glucocorticoids (GC). Has a dual mode of action: as a transcription factor that binds to glucocorticoid response elements (GRE), both for nuclear and mitochondrial DNA, and as a modulator of other transcription factors. Affects inflammatory responses, cellular proliferation and differentiation in target tissues. Involved in chromatin remodeling. Plays a role in rapid mRNA degradation by binding to the 5' UTR of target mRNAs and interacting with PNRC2 in a ligand-dependent manner which recruits the RNA helicase UPF1 and the mRNA-decapping enzyme DCP1A, leading to RNA decay. Could act as a coactivator for STAT5-dependent transcription upon growth hormone (GH) stimulation and could reveal an essential role of hepatic GR in the control of body growth.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Western blot analysis of Glucocorticoid Receptor in rat skeletal muscle lysates using Glucocorticoid Receptor antibody.



Western blot analysis of Glucocorticoid Receptor in K562, 3T3, HeLa lysates using Glucocorticoid Receptor antibody.

产品名称: Glucocorticoid Receptor Rabbit Monoclonal Antibody
产品货号: AMRe02037



Immunocytochemistry analysis of Glucocorticoid Receptor(green) in HeLa using Glucocorticoid Receptor antibody, and DAPI(blue)

注意事项 (Note)

For research use only .