

产品名称: GNB2 Rabbit Monoclonal Antibody
产品货号: AMRe02053



产品概述 (Summary)

产品名称 (Production Name)	GNB2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GNB2
别名 (Alternative Names)	Gnb2; Gnb2l1; RACK1; Transducin beta chain 2
基因 ID (Gene ID)	2783
蛋白 ID (SwissProt ID)	P62879.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa; Observed MW: 32 kDa

研究背景 (Background)

Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane

产品名称: GNB2 Rabbit Monoclonal Antibody
产品货号: AMRe02053




signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of GNB2 in mouse heart lysates using GNB2 antibody.



Western blot analysis of GNB2 in C6, rat Brain, MCF-7 lysates using GNB2 antibody.

注意事项 (Note)

For research use only .