

产品名称: GOLGA1 Rabbit Monoclonal Antibody
产品货号: AMRe02055



产品概述 (Summary)

产品名称 (Production Name)	GOLGA1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GOLGA1
别名 (Alternative Names)	Golgin subfamily A member 1 (Golgin-97)
基因 ID (Gene ID)	2800
蛋白 ID (SwissProt ID)	Q92805.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 88 kDa; Observed MW: 88 kDa

研究背景 (Background)

The Golgi-associated protein golgin A1 (GOLGA1, golgin-97) was first isolated as a Golgi complex autoantigen associated

产品名称: GOLGA1 Rabbit Monoclonal Antibody
产品货号: AMRe02055




with the autoimmune disorder Sjogren's syndrome. The golgin-97 protein contains a carboxy-terminal GRIP domain and is a commonly used trans-Golgi network (TGN) marker.

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of Golgin 97 in HeLa lysates using GOLGA1 antibody.

注意事项 (Note)

For research use only .