

产品名称: HADHA Rabbit Monoclonal Antibody
产品货号: AMRe02069



产品概述 (Summary)

产品名称 (Production Name)	HADHA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HADHA
别名 (Alternative Names)	GBP; ECHA; HADH; LCEH; MTPA; LCHAD; TP-ALPHA
基因 ID (Gene ID)	3030
蛋白 ID (SwissProt ID)	P40939.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 83 kDa; Observed MW: 78 kDa

研究背景 (Background)

Mitochondrial trifunctional enzyme catalyzes the last three of the four reactions of the mitochondrial beta-oxidation

产品名称: HADHA Rabbit Monoclonal Antibody
产品货号: AMRe02069




pathway.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of HADHA in K562, HeLa lysates using HADHA antibody.



Immunocytochemistry analysis of HADHA(green) in HeLa using HADHA antibody, and DAPI(blue)

注意事项 (Note)

For research use only .