

产品名称: HNRNPA1L2 Rabbit Monoclonal Antibody
产品货号: AMRe02078



产品概述 (Summary)

产品名称 (Production Name)	HNRNPA1L2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HNRNPA1L2
别名 (Alternative Names)	hnRNP core protein A1-like 2
基因 ID (Gene ID)	144983.0
蛋白 ID (SwissProt ID)	Q32P51.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 34 kDa

研究背景 (Background)

Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the

产品名称: HNRNPA1L2 Rabbit Monoclonal Antibody
产品货号: AMRe02078




cytoplasm and may modulate splice site selection.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Heterogeneous Nuclear Ribonucleoprotein (A1like) in K562, C6, 3T3, HeLa lysates using HNRNPA1L2 antibody.

注意事项 (Note)

For research use only .