

产品名称: Hexokinase II Rabbit Monoclonal Antibody
产品货号: AMRe02080



产品概述 (Summary)

产品名称 (Production Name)	Hexokinase II Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HK2
别名 (Alternative Names)	HKII; HXK2
基因 ID (Gene ID)	3099
蛋白 ID (SwissProt ID)	P52789.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 102 kDa; Observed MW: 102 kDa

研究背景 (Background)

MiscellaneousIn vertebrates there are four major glucose-phosphorylating isoenzymes, designated hexokinase I, II, III and

产品名称: Hexokinase II Rabbit Monoclonal Antibody
产品货号: AMRe02080




IV (glucokinase).


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of Hexokinase II in rat skeletal muscle lysates using Hexokinase II antibody.



Western blot analysis of Hexokinase II in HeLa, CHO-K1 lysates using Hexokinase II antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Hexokinase II antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Hexokinase II Rabbit Monoclonal Antibody
产品货号: AMRe02080



注意事项 (Note)

For research use only .