

产品名称: HLTF Rabbit Monoclonal Antibody
产品货号: AMRe02096



产品概述 (Summary)

产品名称 (Production Name)	HLTF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HLTF
别名 (Alternative Names)	HIP116; HIP116A; HLTF 1; Hltf; HLTF1; p113; RNF80; SMARCA A3; SMARCA 3; SMARCA3; SNF2L3; ZBU1
基因 ID (Gene ID)	6596
蛋白 ID (SwissProt ID)	Q14527.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 114 kDa; Observed MW: 114 kDa

研究背景 (Background)

产品名称: HLTF Rabbit Monoclonal Antibody
产品货号: AMRe02096




Has both helicase and E3 ubiquitin ligase activities. Possesses intrinsic ATP-dependent nucleosome-remodeling activity; This activity may be required for transcriptional activation or repression of specific target promoters (By similarity).


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of HLTF in HeLa lysates using HLTF antibody.



Immunocytochemistry analysis of HLTF(green) in HeLa using HLTF antibody, and DAPI(blue)

注意事项 (Note)

For research use only .