

产品名称: HNF 4 alpha Rabbit Monoclonal Antibody
产品货号: AMRe02102



产品概述 (Summary)

产品名称 (Production Name)	HNF 4 alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP,ChIP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Hnf4A
别名 (Alternative Names)	Hnf4; HNF-4; MODY1; Nr2a1; Tcf14; TCF-14; Hnf4alpha
基因 ID (Gene ID)	15378.0
蛋白 ID (SwissProt ID)	P49698.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50,ChIP 1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 53 kDa; Observed MW: 53 kDa

研究背景 (Background)

Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin,

产品名称: HNF 4 alpha Rabbit Monoclonal Antibody
产品货号: AMRe02102




apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and inTine.MiscellaneousBinds fatty acids.


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of HNF4alpha in mouse liver lysates using HNF 4 alpha antibody.



Western blot analysis of HNF4alpha in K562, C6, HeLa lysates using HNF4alpha antibody.

注意事项 (Note)

For research use only .