

产品名称: HOXB4 Rabbit Monoclonal Antibody
产品货号: AMRe02112



产品概述 (Summary)

产品名称 (Production Name)	HOXB4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HOXB4
别名 (Alternative Names)	HOXB4; HOX2F; Homeobox protein Hox-B4; Homeobox protein Hox-2.6; Homeobox protein Hox-2F
基因 ID (Gene ID)	3214
蛋白 ID (SwissProt ID)	P17483.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 34 kDa

研究背景 (Background)

产品名称: HOXB4 Rabbit Monoclonal Antibody
产品货号: AMRe02112




Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Immunocytochemistry analysis of HOXB4 (green) in K562 using HOXB4 antibody, and DAPI (blue).



Western blot analysis of HOXB4 in K562 lysates using HOXB4 antibody.

注意事项 (Note)

For research use only .