

产品名称: HRH3 Rabbit Monoclonal Antibody
产品货号: AMRe02118



产品概述 (Summary)

产品名称 (Production Name)	HRH3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HRH3
别名 (Alternative Names)	HRH3; GPCR97; Histamine H3 receptor; H3R; HH3R; G-protein coupled receptor 97
基因 ID (Gene ID)	11255
蛋白 ID (SwissProt ID)	Q9Y5N1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 49 kDa; Observed MW: 49 kDa

研究背景 (Background)

产品名称: HRH3 Rabbit Monoclonal Antibody
产品货号: AMRe02118




The H3 subclass of histamine receptors could mediate the histamine signals in CNS and peripheral nervous system. Signals through the inhibition of adenylate cyclase and displays high constitutive activity (spontaneous activity in the absence of agonist).

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of HRH3/Histamine 3 Receptor in K562, C6, 3T3 lysates using HRH3 antibody.

注意事项 (Note)

For research use only .