

产品名称: Hsp90 alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe02127



产品概述 (Summary)

产品名称 (Production Name)	Hsp90 alpha/beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant Rabbit Monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC-F,IHC-P,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal Antibody
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
纯化方式 (Purification)	Affinity Purified

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSP90AA1/HSP90AB1
别名 (Alternative Names)	HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38; HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
基因 ID (Gene ID)	3320/3326
蛋白 ID (SwissProt ID)	P07900/P08238

产品应用 (Application)

稀释比 (Dilution Ratio)	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20
蛋白分子量 (Molecular Weight)	Calculated MW: 85 kDa; Observed MW: 90 kDa

产品名称: Hsp90 alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe02127




研究背景 (Background)

Molecular chaperone that promotes the maturation, structural maintenance and proper regulation of specific target proteins involved for instance in cell cycle control and signal transduction. Undergoes a functional cycle that is linked to its ATPase activity. This cycle probably induces conformational changes in the client proteins, thereby causing their activation.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Western blot analysis of Hsp90 alpha/beta in rat Brain, C6, 3T3, HeLa lysates using Hsp90 alpha/beta antibody.



Immunocytochemistry analysis of Hsp90 alpha/beta(green) in HeLa using Hsp90 alpha/beta antibody, and DAPI(blue)

产品名称: Hsp90 alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe02127



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Hsp90 alpha/beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .