

产品名称: HuC/HuD Rabbit Monoclonal Antibody
产品货号: AMRe02130



产品概述 (Summary)

产品名称 (Production Name)	HuC/HuD Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Elavl3/Elavl4
别名 (Alternative Names)	Huc; mHuC; PLE21; 2600009P04Rik
基因 ID (Gene ID)	15571/1996
蛋白 ID (SwissProt ID)	Q60900/P26378.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa; Observed MW: 40 kDa

研究背景 (Background)

Binds to AU-rich sequences (AREs) of target mRNAs, including VEGF mRNA. May also bind poly-A tracts via RRM 3. May be

产品名称: HuC/HuD Rabbit Monoclonal Antibody
产品货号: AMRe02130




involved in neuronal differentiation and maintenance.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of HuD/HuC in rat Brain lysates using HuC/HuD antibody.

注意事项 (Note)

For research use only .