

产品名称: ICOS Ligand Rabbit Monoclonal Antibody
产品货号: AMRe02132



产品概述 (Summary)

产品名称 (Production Name)	ICOS Ligand Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ICOSLG
别名 (Alternative Names)	inducible T-cell co-stimulator ligand; B7H2; GL50; B7-H2; B7RP1; CD275; ICOSL; LICOS; B7RP-1; ICOS-L
基因 ID (Gene ID)	23308
蛋白 ID (SwissProt ID)	O75144.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 75 kDa

研究背景 (Background)

产品名称: ICOS Ligand Rabbit Monoclonal Antibody
产品货号: AMRe02132




Ligand for the T-cell-specific cell surface receptor ICOS. Acts as a costimulatory signal for T-cell proliferation and cytokine secretion; induces also B-cell proliferation and differentiation into plasma cells. Could play an important role in mediating local tissue responses to inflammatory conditions, as well as in modulating the secondary immune response by co-stimulating memory T-cell function .

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of ICOS Ligand in K562 lysates using ICOS Ligand antibody.

注意事项 (Note)

For research use only .