

产品名称: INMT Rabbit Monoclonal Antibody
产品货号: AMRe02158



产品概述 (Summary)

产品名称 (Production Name)	INMT Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	INMT
别名 (Alternative Names)	TEMT
基因 ID (Gene ID)	11185
蛋白 ID (SwissProt ID)	O95050.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 29 kDa

研究背景 (Background)

Functions as thioether S-methyltransferase and is active with a variety of thioethers and the corresponding selenium and

产品名称: INMT Rabbit Monoclonal Antibody
产品货号: AMRe02158




tellurium compounds, including 3-methylthiopropionaldehyde, dimethyl selenide, dimethyl telluride, 2-methylthioethylamine, 2-methylthioethanol, methyl-n-propyl sulfide and diethyl sulfide. Plays an important role in the detoxification of selenium compounds . Catalyzes the N-methylation of tryptamine and structurally related compounds.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of INMT in 3T3 lysates using INMT antibody.

注意事项 (Note)

For research use only .