

产品名称: IQGAP2 Rabbit Monoclonal Antibody
产品货号: AMRe02168



产品概述 (Summary)

产品名称 (Production Name)	IQGAP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IQGAP2
别名 (Alternative Names)	Ras GTPase-activating-like protein IQGAP2
基因 ID (Gene ID)	10788
蛋白 ID (SwissProt ID)	Q13576.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 181 kDa; Observed MW: 181 kDa

研究背景 (Background)

Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity. Associates with calmodulin.


产品名称: IQGAP2 Rabbit Monoclonal Antibody
产品货号: AMRe02168



研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of IQGAP2 in Ramos lysates using IQGAP2 antibody.

注意事项 (Note)

For research use only .