

产品名称: IRAK4 Rabbit Monoclonal Antibody
产品货号: AMRe02169



产品概述 (Summary)

产品名称 (Production Name)	IRAK4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IRAK4
别名 (Alternative Names)	IL-1 receptor-associated kinase 4; Interleukin 1 receptor associated kinase 4 mutant form 1; IPD1; IRAK4; LOC 51135; NY REN 64; REN64
基因 ID (Gene ID)	51135
蛋白 ID (SwissProt ID)	Q9NWZ3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 52 kDa

研究背景 (Background)

产品名称: IRAK4 Rabbit Monoclonal Antibody
产品货号: AMRe02169




Required for the efficient recruitment of IRAK1 to the IL-1 receptor complex following IL-1 engagement, triggering intracellular signaling cascades leading to transcriptional up-regulation and mRNA stabilization. Phosphorylates IRAK1.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Western blot analysis of IRAK4 in Jurkat, HeLa lysates using IRAK4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IRAK4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of IRAK4 (green) in HeLa using IRAK4 antibody, and DAPI (blue)

产品名称: IRAK4 Rabbit Monoclonal Antibody
产品货号: AMRe02169



注意事项 (Note)

For research use only .