

产品名称: JunB Rabbit Monoclonal Antibody
产品货号: AMRe02180



产品概述 (Summary)

产品名称 (Production Name)	JunB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	JUNB
别名 (Alternative Names)	Transcription factor jun-B
基因 ID (Gene ID)	3726
蛋白 ID (SwissProt ID)	P17275.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 36 kDa; Observed MW: 42,43 kDa

研究背景 (Background)

JunB is a basic region, leucine zipper (bZIP) transcription factor belonging to the Jun family that includes c-Jun and JunD.

产品名称: JunB Rabbit Monoclonal Antibody
产品货号: AMRe02180




Jun family members homodimerize or heterodimerize with Fos and ATF proteins to form a functional transcription factor AP-1 (activator protein, whose activity is regulated by a variety of physiological and pathological stimuli such as growth factors, infections, and stress signals).

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of JunB in K562, C6, 3T3, HeLa lysates using JunB antibody.

注意事项 (Note)

For research use only .